

产品名称: RBM3 (18E7) Rabbit Monoclonal Antibody
产品货号: AMRe16955



产品概述 (Summary)

产品名称 (Production Name)	RBM3 (18E7) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBM3
别名 (Alternative Names)	Rbm3; RNPL;
基因 ID (Gene ID)	5935.0
蛋白 ID (SwissProt ID)	P98179.A synthetic peptide of human RBM3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:20-1:50, FC 1:20-1:50
蛋白分子量 (Molecular Weight)	17kDa

研究背景 (Background)

Cold-inducible mRNA binding protein that enhances global protein synthesis at both physiological and mild hypothermic

产品名称: RBM3 (18E7) Rabbit Monoclonal Antibody
产品货号: AMRe16955




temperatures. Reduces the relative abundance of microRNAs, when overexpressed. Enhances phosphorylation of translation initiation factors and active polysome formation. Cold-inducible mRNA binding protein that enhances global protein synthesis at both physiological and mild hypothermic temperatures. Reduces the relative abundance of microRNAs, when overexpressed. Enhances phosphorylation of translation initiation factors and active polysome formation (By similarity).

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot detection of RBM3 in A549,HL-60,U251 cell lysates using RBM3 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .