

产品名称: RCHY1 (11B14) Rabbit Monoclonal Antibody
产品货号: AMRe16981



产品概述 (Summary)

产品名称 (Production Name)	RCHY1 (11B14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RCHY1
别名 (Alternative Names)	ARNIP; CHIMP; hARNIP; hPirh2; PIRH2E; PIRH2F; RCHY1; RNF199; ZCHY; ZNF363;
基因 ID (Gene ID)	25898.0
蛋白 ID (SwissProt ID)	Q96PM5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:100-1:200
蛋白分子量 (Molecular Weight)	30kDa

产品名称: RCHY1 (11B14) Rabbit Monoclonal Antibody
产品货号: AMRe16981




研究背景 (Background)

Mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including p53/TP53, HDAC1 and CDKN1B. Preferentially acts on tetrameric p53/TP53. Mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including p53/TP53, P73, HDAC1 and CDKN1B. Preferentially acts on tetrameric p53/TP53. Monoubiquitinates the translesion DNA polymerase POLH. Contributes to the regulation of the cell cycle progression. Increases AR transcription factor activity.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of RCHY1 expression in 293 cell lysate.

注意事项 (Note)

For research use only .