

产品名称: REA (13I14) Rabbit Monoclonal Antibody
产品货号: AMRe16988



产品概述 (Summary)

产品名称 (Production Name)	REA (13I14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PHB2
别名 (Alternative Names)	BAP; Bap37; BCAP 37; D prohibitin; p22; Phb2; PNAS 141; Prohibitin 2;
基因 ID (Gene ID)	11331.0
蛋白 ID (SwissProt ID)	Q99623.A synthetic peptide of human REA

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:1000-1:5000,ICC/IF 1:20-1:50,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	33kDa

研究背景 (Background)

Acts as a mediator of transcriptional repression by nuclear hormone receptors via recruitment of histone deacetylases (By

产品名称: REA (13I14) Rabbit Monoclonal Antibody
产品货号: AMRe16988




similarity). Functions as an estrogen receptor (ER)-selective coregulator that potentiates the inhibitory activities of antiestrogens and represses the activity of estrogens. Protein with pleiotropic attributes mediated in a cell- compartment- and tissue-specific manner, which include the plasma membrane-associated cell signaling functions, mitochondrial chaperone, and transcriptional co-regulator of transcription factors and sex steroid hormones in the nucleus.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of REA in K562,Rat Brain,C6,3T3,HeLa cell lysates using REA antibody(1:1000 diluted).

注意事项 (Note)

For research use only .