

产品名称: RGS6 (3A18) Rabbit Monoclonal Antibody
产品货号: AMRe17100



产品概述 (Summary)

产品名称 (Production Name)	RGS6 (3A18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RGS6
别名 (Alternative Names)	GAP; RGS6; S914;
基因 ID (Gene ID)	9628.0
蛋白 ID (SwissProt ID)	P49758.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	54kDa

研究背景 (Background)


产品名称: RGS6 (3A18) Rabbit Monoclonal Antibody
产品货号: AMRe17100



Inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits thereby driving them into their inactive GDP-bound form. Activity on G(o)-alpha is specifically enhanced by the RGS6/Gbeta5 dimer. Regulates G protein-coupled receptor signaling cascades. Inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits, thereby driving them into their inactive GDP-bound form. The RGS6/GNB5 dimer enhances GNAO1 GTPase activity (PubMed: [10521509](http://www.uniprot.org/citations/10521509)).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RGS6 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .