

产品名称: RHEB (18Y18) Rabbit Monoclonal Antibody
产品货号: AMRe17106



产品概述 (Summary)

产品名称 (Production Name)	RHEB (18Y18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RHEB
别名 (Alternative Names)	RHEB 2; Rheb; RHEB2;
基因 ID (Gene ID)	6009.0
蛋白 ID (SwissProt ID)	Q15382.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:100-1:200,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	20kDa

研究背景 (Background)


产品名称: RHEB (18Y18) Rabbit Monoclonal Antibody
产品货号: AMRe17106



Stimulates the phosphorylation of S6K1 and EIF4EBP1 through activation of mTORC1 signaling. Activates the protein kinase activity of mTORC1. Has low intrinsic GTPase activity. Activates the protein kinase activity of mTORC1, and thereby plays a role in the regulation of apoptosis. Stimulates the phosphorylation of S6K1 and EIF4EBP1 through activation of mTORC1 signaling. Has low intrinsic GTPase activity.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RHEB expression in Raji cell lysate.

注意事项 (Note)

For research use only .