

产品名称: RPS20 (7B14) Rabbit Monoclonal Antibody
产品货号: AMRe17372



产品概述 (Summary)

产品名称 (Production Name)	RPS20 (7B14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPS20
别名 (Alternative Names)	40S ribosomal protein S20; RS20; Ribosomal protein S20; RPS 20;
基因 ID (Gene ID)	6224.0
蛋白 ID (SwissProt ID)	P60866.A synthetic peptide of human RPS20

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:100,ICC/IF 1:20-1:50,FC 1:20-1:50
蛋白分子量 (Molecular Weight)	13kDa

研究背景 (Background)

RPS20 is a a ribosomal protein that is a component of the 40S subunit. Ribosomes, the organelles that catalyze protein

产品名称: RPS20 (7B14) Rabbit Monoclonal Antibody
产品货号: AMRe17372




synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. RPS20 belongs to the S10P family of ribosomal proteins and is located in the cytoplasm.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot detection of RPS20 in HeLa,CHO-K1,C6,Ramos using RPS20 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .