

产品名称: RSK4 (7016) Rabbit Monoclonal Antibody
产品货号: AMRe17410



产品概述 (Summary)

产品名称 (Production Name)	RSK4 (7016) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPS6KA6
别名 (Alternative Names)	p90RSK6; pp90RSK4; Ribosomal S6 kinase 4; rps6ka6; RSK 4;
基因 ID (Gene ID)	27330.0
蛋白 ID (SwissProt ID)	Q9UK32.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:100-1:1000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	84kDa

研究背景 (Background)

产品名称: RSK4 (7016) Rabbit Monoclonal Antibody
产品货号: AMRe17410




Serine/threonine kinase that may play a role in mediating the growth-factor and stress induced activation of the transcription factor CREB. Constitutively active serine/threonine-protein kinase that exhibits growth-factor-independent kinase activity and that may participate in p53/TP53-dependent cell growth arrest signaling and play an inhibitory role during embryogenesis.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of RSK4 expression in SH-SY5Y cell lysate.

注意事项 (Note)

For research use only .