

产品名称: Salivary alpha amylase (5D14) Rabbit Monoclonal Antibody
产品货号: AMRe17579



产品概述 (Summary)

产品名称 (Production Name)	Salivary alpha amylase (5D14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	
别名 (Alternative Names)	Alpha amylase; AMY1; AMY1C; AMY2; AMY2B; AMY3; HXA;
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P04745.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	58kDa

产品名称: Salivary alpha amylase (5D14) Rabbit Monoclonal Antibody
产品货号: AMRe17579




研究背景 (Background)

Belongs to the glycosyl hydrolase 13 family.

研究领域 (Research Area)

Starch and sucrose metabolism;

图片 (Image Data)



Western blot analysis of Salivary alpha amylase expression in Human pancreas lysate.

注意事项 (Note)

For research use only .