

产品名称: SEC61A (18P13) Rabbit Monoclonal Antibody
产品货号: AMRe17694



产品概述 (Summary)

产品名称 (Production Name)	SEC61A (18P13) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	
别名 (Alternative Names)	fb62c11; HSEC61; SEC 61; SEC 61A; Sec61 alpha 1; SEC61; SEC61A1; SEC61A2;
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P61619/Q9H9S3.Recombinant protein of human SEC61A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:20-1:50,ICC/IF 1:20-1:50,FC 1:20-1:50,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	52kDa

研究背景 (Background)

Plays a crucial role in the insertion of secretory and membrane polypeptides into the ER. Required for assembly of


产品名称: SEC61A (18P13) Rabbit Monoclonal Antibody
产品货号: AMRe17694



membrane and secretory proteins. Tightly associated with membrane-bound ribosomes, either directly or through adapter proteins.

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of SEC61A in K562,C6,3T3 cell lysates using SEC61A antibody(1:1000 diluted).

注意事项 (Note)

For research use only .