

**产品名称: SERCA2 (13M14) Rabbit Monoclonal Antibody**

**产品货号: AMRe17748**

### 产品概述 (Summary)

|                   |                                                                                                                             |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------|
| 描述(Description)   | Recombinant rabbit monoclonal antibody                                                                                      |
| 宿主(Host)          | Rabbit                                                                                                                      |
| 应用(Application)   | WB,IHC,ICC/IF,FC                                                                                                            |
| 种属反应(Reactivity)  | Human,Mouse,Rat                                                                                                             |
| 偶联物(Conjugation)  | Unconjugated                                                                                                                |
| 修饰(Modification)  | Unmodified                                                                                                                  |
| 同种型(Isotype)      | IgG                                                                                                                         |
| 克隆(Clonality)     | Monoclonal                                                                                                                  |
| 剂型(Form)          | Liquid                                                                                                                      |
| 保存条件(Storage)     | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                    |
| 储存溶液(Buffer)      | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein. |
| 纯化 (Purification) | Affinity purification                                                                                                       |

### 产品应用(Application)

|                       |                                                                 |
|-----------------------|-----------------------------------------------------------------|
| 稀释比(Dilution Ratio)   | WB 1:2000-1:10000,IHC 1:20-1:100,ICC/IF 1:20-1:50,FC 1:20-1:100 |
| 分子量(Molecular Weight) | 115kDa                                                          |

### 抗原信息(Antigen Information)

|                       |                                            |
|-----------------------|--------------------------------------------|
| 基因名(Gene Name)        | ATP2A2                                     |
| 别名(Alternative Names) | SERCA2; Calcium pump 2; ATP2B;             |
| 基因 ID(Gene ID)        | 488.0                                      |
| SwissProt ID          | P16615                                     |
| 免疫原(Immunogen)        | A synthetic peptide of human SERCA2 ATPase |

### 研究背景 (Background)

This magnesium-dependent enzyme catalyzes the hydrolysis of ATP coupled with the translocation of calcium from the cytosol to the sarcoplasmic reticulum lumen. Isoform 2 is involved in the regulation of the contraction/relaxation cycle. This magnesium-dependent enzyme catalyzes the hydrolysis of ATP coupled with the translocation of calcium from the cytosol to

the sarcoplasmic reticulum lumen (PubMed:16402920). Involved in autophagy in response to starvation. Upon interaction with VMP1 and activation, controls ER-isolation membrane contacts for autophagosome formation (PubMed:28890335). Also modulates ER contacts with lipid droplets, mitochondria and endosomes (PubMed:28890335).


## 研究领域 (Research Area)

Neuroscience

## 注意事项 (Note)

For research use only.

## 图片 (Image Data)



Western blot detection of SERCA2 ATPase in K562,Hela cell lysates using SERCA2 ATPase antibody(1:1000 diluted).

## EnkiLife 优势产品推荐

| 产品类别                        | 产品货号     | 产品名称                                               |
|-----------------------------|----------|----------------------------------------------------|
| WB 解决方案                     | RA10020  | 2h 极速 WB 即用型全流程试剂盒                                 |
|                             | RA10021  | 4h 快速 WB 即用型全流程试剂盒                                 |
|                             | RA10042  | 5h 畅享版 WB 全流程试剂盒                                   |
|                             | RA10037  | 校准级彩色预染蛋白 Marker (8-180kDa)                        |
|                             | RA10038  | 校准级彩色预染蛋白 Marker (10-250kDa)                       |
|                             | RA10039  | 校准级高分子彩色预染蛋白 Marker (25-400kDa)                    |
| TSA 多重荧光染色试剂盒               | RA10008  | TSA 双标三色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10009  | TSA 三标四色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10010  | TSA 四标五色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10011  | TSA 五标六色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10012  | TSA 六标七色多重荧光染色试剂盒 (mIHC)                           |
| IHC 检测试剂盒                   | RA10006  | HRP Anti-Mouse/Rabbit IHC Detection System         |
|                             | RA10007  | Polymer-HRP Anti-Mouse/Rabbit IHC Detection System |
| 抗体标记试剂盒                     | RE80004p | 辣根过氧化物酶(HRP)抗体标记试剂盒                                |
|                             | RE80002q | Sulfo-NHS-生物素标记试剂盒                                 |
|                             | RE80007p | Cy3 荧光素标记试剂盒                                       |
|                             | RE80011p | Fluor488 荧光素标记试剂盒                                  |
|                             | RE80017p | Fluor750 荧光素标记试剂盒                                  |
|                             | RE80005p | 藻红蛋白(R-PE) 抗体快速标记试剂盒                               |
|                             | RE80040  | PE-Cy7 串联染料抗体快速标记试剂盒                               |
| 稳转细胞系构建服务<br>(免费赠送全膜 WB 验证) | TS-0001  | 过表达稳转细胞系构建                                         |
|                             | TS-0002  | 敲低稳转细胞系构建                                          |
|                             | TS-0003  | 敲除细胞系构建                                            |

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>