

产品名称: SERPING1 (12Y5) Rabbit Monoclonal Antibody
产品货号: AMRe17768



产品概述 (Summary)

产品名称 (Production Name)	SERPING1 (12Y5) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SERPING1
别名 (Alternative Names)	C1IN; C1Inh; C1NH; HAE1; HAE2;
基因 ID (Gene ID)	710.0
蛋白 ID (SwissProt ID)	P05155.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	55kDa

研究背景 (Background)

产品名称: SERPING1 (12Y5) Rabbit Monoclonal Antibody
产品货号: AMRe17768




Activation of the C1 complex is under control of the C1-inhibitor. It forms a proteolytically inactive stoichiometric complex with the C1r or C1s proteases. Activation of the C1 complex is under control of the C1- inhibitor. It forms a proteolytically inactive stoichiometric complex with the C1r or C1s proteases. May play a potentially crucial role in regulating important physiological pathways including complement activation, blood coagulation, fibrinolysis and the generation of kinins. Very efficient inhibitor of FXIIa. Inhibits chymotrypsin and kallikrein.

研究领域 (Research Area)

Complement and coagulation cascades;

图片 (Image Data)



Western blot analysis of SERPING1 expression in Human plasma lysate.

注意事项 (Note)

For research use only .