

产品名称: SLC12A1 / NKCC2 (7C8) Rabbit Monoclonal Antibody
产品货号: AMRe17952



产品概述 (Summary)

产品名称 (Production Name)	SLC12A1 / NKCC2 (7C8) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC12A1
别名 (Alternative Names)	BSC1; NKCC2; Slc12a1;
基因 ID (Gene ID)	6557.0
蛋白 ID (SwissProt ID)	Q13621.Recombinant protein of human NKCC2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	122kDa

研究背景 (Background)


产品名称: SLC12A1 / NKCC2 (7C8) Rabbit Monoclonal Antibody
产品货号: AMRe17952

May play an important role in cell growth and/or transcription. Renal sodium, potassium and chloride ion cotransporter that mediates the transepithelial NaCl reabsorption in the thick ascending limb and plays an essential role in the urinary concentration and volume regulation. Electrically silent transporter system.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot detection of SLC12A1/NKCC2 in HeLa cell lysates using SLC12A1/NKCC2 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .