

产品名称: SNAP29 (10N2) Rabbit Monoclonal Antibody
产品货号: AMRe18048



产品概述 (Summary)

产品名称 (Production Name)	SNAP29 (10N2) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SNAP29
别名 (Alternative Names)	CEDNIK; SNAP29;
基因 ID (Gene ID)	9342.0
蛋白 ID (SwissProt ID)	O95721.A synthetic peptide of human SNAP29

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:20-1:100, FC 1:20-1:50, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	29kDa

研究背景 (Background)

Involved in multiple membrane trafficking steps. SNAREs, soluble N-ethylmaleimide-sensitive factor-attachment protein

产品名称: SNAP29 (10N2) Rabbit Monoclonal Antibody
产品货号: AMRe18048




receptors, are essential proteins for fusion of cellular membranes. SNAREs localized on opposing membranes assemble to form a trans-SNARE complex, an extended, parallel four alpha-helical bundle that drives membrane fusion. SNAP29 is a SNARE involved in autophagy through the direct control of autophagosome membrane fusion with the lysosome membrane. Plays also a role in ciliogenesis by regulating membrane fusions.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of extracts from HeLa cells, Mouse brain and Rat brain tissues using RM6766 at 1:1000.

注意事项 (Note)

For research use only .