

产品名称: STAT2 (5I18) Rabbit Monoclonal Antibody
产品货号: AMRe18349



产品概述 (Summary)

产品名称 (Production Name)	STAT2 (5I18) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STAT2
别名 (Alternative Names)	
基因 ID (Gene ID)	6773.0
蛋白 ID (SwissProt ID)	P52630.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200,ICC/IF 1:20-1:50,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	98 kDa

研究背景 (Background)


产品名称: STAT2 (5I18) Rabbit Monoclonal Antibody
产品货号: AMRe18349




研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of STAT2 in K562 lysates using STAT2 antibody.




Immunohistochemistry analysis of paraffin-embedded Human lung cancer using STAT2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of STAT2 in A431 using STAT2 antibody.

产品名称: STAT2 (5I18) Rabbit Monoclonal Antibody
产品货号: AMRe18349



Western blot analysis of STAT2 in HeLa , Rat brain lysates using STAT2 antibody.

注意事项 (Note)

For research use only .