

产品名称: Staufen (11R2) Rabbit Monoclonal Antibody  
产品货号: AMRe18374



## 产品概述 (Summary)

|                        |                                           |
|------------------------|-------------------------------------------|
| 产品名称 (Production Name) | Staufen (11R2) Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody    |
| 宿主 (Host)              | Rabbit                                    |
| 应用 (Application)       | WB,IHC,ICC/IF,FC                          |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                           |

## 产品性能 (Performance)

|                     |                                                                                                                             |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                |
| 修饰 (Modification)   | Unmodified                                                                                                                  |
| 同种型 (Isotype)       | IgG                                                                                                                         |
| 克隆 (Clonality)      | Monoclonal                                                                                                                  |
| 形式 (Form)           | Liquid                                                                                                                      |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                    |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein. |
| 纯化方式 (Purification) | Affinity purification                                                                                                       |

## 免疫原信息 (Immunogen)

|                        |                                                                                        |
|------------------------|----------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | STAU1<br>Double-stranded RNA-binding protein Staufen homolog 1; PPP1R150; STAU;        |
| 别名 (Alternative Names) | staufen; Staufen RNA binding protein (Drosophila); Staufen, Drosophila, homolog of, 1; |
| 基因 ID (Gene ID)        | 6780.0                                                                                 |
| 蛋白 ID (SwissProt ID)   | O95793.A synthetic peptide of human Staufen                                            |

## 产品应用 (Application)

|                          |                                                                   |
|--------------------------|-------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:100-1:200,FC 1:200-1:500 |
| 蛋白分子量 (Molecular Weight) | 63kDa                                                             |

产品名称: Staufen (11R2) Rabbit Monoclonal Antibody  
产品货号: AMRe18374



## 研究背景 (Background)


Binds double-stranded RNA (regardless of the sequence) and tubulin. May play a role in specific positioning of mRNAs at given sites in the cell by cross-linking cytoskeletal and RNA components, and in stimulating their translation at the site.

Binds double-stranded RNA (regardless of the sequence) and tubulin. May play a role in specific positioning of mRNAs at given sites in the cell by cross-linking cytoskeletal and RNA components, and in stimulating their translation at the site.

## 研究领域 (Research Area)

Signal Transduction

## 图片 (Image Data)



Western blot detection of Staufen in K562,C6,3T3 cell lysates using Staufen antibody(1:1000 diluted).

## 注意事项 (Note)

For research use only .