

产品名称: Syntaxin 16 (10E13) Rabbit Monoclonal Antibody  
产品货号: AMRe18512



## 产品概述 (Summary)

|                        |                                                |
|------------------------|------------------------------------------------|
| 产品名称 (Production Name) | Syntaxin 16 (10E13) Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody         |
| 宿主 (Host)              | Rabbit                                         |
| 应用 (Application)       | WB, ICC/IF                                     |
| 种属反应性 (Reactivity)     | Human, Mouse, Rat                              |

## 产品性能 (Performance)

|                     |                                                                                                                             |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                |
| 修饰 (Modification)   | Unmodified                                                                                                                  |
| 同种型 (Isotype)       | IgG                                                                                                                         |
| 克隆 (Clonality)      | Monoclonal                                                                                                                  |
| 形式 (Form)           | Liquid                                                                                                                      |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                    |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein. |
| 纯化方式 (Purification) | Affinity purification                                                                                                       |

## 免疫原信息 (Immunogen)

|                        |                                                 |
|------------------------|-------------------------------------------------|
| 基因名 (Gene Name)        | STX16                                           |
| 别名 (Alternative Names) | hsyn16; Stx16; Syn16;                           |
| 基因 ID (Gene ID)        | 8675.0                                          |
| 蛋白 ID (SwissProt ID)   | O14662.A synthetic peptide of human Syntaxin 16 |

## 产品应用 (Application)

|                          |                                      |
|--------------------------|--------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000, ICC/IF 1:100-1:200 |
| 蛋白分子量 (Molecular Weight) | 37kDa                                |

## 研究背景 (Background)

SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network. SNARE involved in vesicular

产品名称: Syntaxin 16 (10E13) Rabbit Monoclonal Antibody  
产品货号: AMRe18512




transport from the late endosomes to the trans-Golgi network.

## 研究领域 (Research Area)

Tags & Cell Markers

## 图片 (Image Data)



Western blot detection of Syntaxin 16 in K562,Rat Brain lysates using Syntaxin 16 antibody(1:1000 diluted).

## 注意事项 (Note)

For research use only .