

产品名称: TAF3 (4T14) Rabbit Monoclonal Antibody
产品货号: AMRe18619



产品概述 (Summary)

产品名称 (Production Name)	TAF3 (4T14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAF3
别名 (Alternative Names)	TAF140; TAF3; TAFII140; TAFII155 protein;
基因 ID (Gene ID)	83860.0
蛋白 ID (SwissProt ID)	Q5VWG9.Recombinant protein of human TAF3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:20-1:100
蛋白分子量 (Molecular Weight)	104kDa

研究背景 (Background)

Transcription factor TFIID is one of the general factors required for accurate and regulated initiation by RNA polymerase II.

产品名称: TAF3 (4T14) Rabbit Monoclonal Antibody
产品货号: AMRe18619




TFIID is a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. Required in complex with TBPL2 for the differentiation of myoblasts into myocytes. The complex replaces TFIID at specific promoters at an early stage in the differentiation process. Transcription factor TFIID is one of the general factors required for accurate and regulated initiation by RNA polymerase II. TFIID is a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. Required in complex with TBPL2 for the differentiation of myoblasts into myocytes. The complex replaces TFIID at specific promoters at an early stage in the differentiation process.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot detection of TAF3 in Jurkat cell lysates using TAF3 antibody(1:1000 diluted).

注意事项 (Note)

For research use only .