

产品名称: Trefoil Factor 3 (11U9) Rabbit Monoclonal Antibody
产品货号: AMRe19232



产品概述 (Summary)

产品名称 (Production Name)	Trefoil Factor 3 (11U9) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP,IF-P
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TFF3
别名 (Alternative Names)	hITF; hP1.B; Intestinal trefoil factor; ITF; mITF; P1B; Polypeptide P1.B; TFF3; TFI; TREFOIL; Trefoil factor 3;
基因 ID (Gene ID)	7033.0
蛋白 ID (SwissProt ID)	Q07654.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:2000,IP 1:20-1:50,IF-P 1:200-1:2000
蛋白分子量 (Molecular Weight)	10kDa

产品名称: Trefoil Factor 3 (11U9) Rabbit Monoclonal Antibody
产品货号: AMRe19232




研究背景 (Background)

Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen). Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Trefoil Factor 3 expression in human small intestine lysate.

注意事项 (Note)

For research use only .