

产品名称: TSFM (9P10) Rabbit Monoclonal Antibody
产品货号: AMRe19356



产品概述 (Summary)

产品名称 (Production Name)	TSFM (9P10) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TSFM
别名 (Alternative Names)	COXPD3; EFTS; EFTSMT; TSFM; TSMT;
基因 ID (Gene ID)	10102.0
蛋白 ID (SwissProt ID)	P43897.A synthetic peptide of human TSFM

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000, ICC/IF 1:50-1:100
蛋白分子量 (Molecular Weight)	35kDa

研究背景 (Background)

Associates with the EF-Tu.GDP complex and induces the exchange of GDP to GTP. It remains bound to the aminoacyl-

产品名称: TSFM (9P10) Rabbit Monoclonal Antibody
产品货号: AMRe19356




tRNA.EF-Tu.GTP complex up to the GTP hydrolysis stage on the ribosome. Associates with the EF-Tu.GDP complex and induces the exchange of GDP to GTP. It remains bound to the aminoacyl-tRNA.EF- Tu.GTP complex up to the GTP hydrolysis stage on the ribosome.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot detection of TSFM in HeLa,A549,HL-60,U2OS cell lysates using TSFM antibody(1:1000 diluted).

注意事项 (Note)

For research use only .