

产品名称: ULK1 (3G10) Rabbit Monoclonal Antibody
产品货号: AMRe19620



产品概述 (Summary)

产品名称 (Production Name)	ULK1 (3G10) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ULK1
别名 (Alternative Names)	
基因 ID (Gene ID)	8408.0
蛋白 ID (SwissProt ID)	O75385.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000
蛋白分子量 (Molecular Weight)	113 kDa

研究背景 (Background)


产品名称: ULK1 (3G10) Rabbit Monoclonal Antibody
产品货号: AMRe19620



研究领域 (Research Area)

Signal Transduction; Protein Phosphorylation; Ser / Thr Kinases; Neuroscience; Neurology process; Neurogenesis; Cardiovascular; Heart; Autophagy; Metabolism; Pathways and Processes; Mitochondrial Metabolism; Mitophagy fission and fusion

图片 (Image Data)



Western blot analysis of ULK1 in HEK293 , PC-12 lysates using ULK1 antibody.

注意事项 (Note)

For research use only .