

产品名称: VEGFD (7C13) Rabbit Monoclonal Antibody
产品货号: AMRe19774



产品概述 (Summary)

产品名称 (Production Name)	VEGFD (7C13) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VEGFD
别名 (Alternative Names)	FIGF; VEGFD;
基因 ID (Gene ID)	2277.0
蛋白 ID (SwissProt ID)	O43915.A synthetic peptide of human VEGFD

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,ICC/IF 1:20-1:50,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	40kDa

研究背景 (Background)

Growth factor active in angiogenesis, lymphangiogenesis and endothelial cell growth, stimulating their proliferation and

产品名称: VEGFD (7C13) Rabbit Monoclonal Antibody
产品货号: AMRe19774




migration and also has effects on the permeability of blood vessels. Growth factor active in angiogenesis, lymphangiogenesis and endothelial cell growth, stimulating their proliferation and migration and also has effects on the permeability of blood vessels. May function in the formation of the venous and lymphatic vascular systems during embryogenesis, and also in the maintenance of differentiated lymphatic endothelium in adults. Binds and activates VEGFR-2 (KDR/FLK1) and VEGFR-3 (FLT4) receptors.

研究领域 (Research Area)

Cytokine-cytokine receptor interaction;mTOR;Focal adhesion;Pathways in cancer;Renal cell carcinoma;Pancreatic cancer;Bladder cancer;

图片 (Image Data)



Western blot detection of VEGFD in Rat Brain,C6,Hela cell lysates using VEGFD antibody(1:1000 diluted).

注意事项 (Note)

For research use only .