

产品名称: Vinculin (19R11) Rabbit Monoclonal Antibody
产品货号: AMRe19800



产品概述 (Summary)

产品名称 (Production Name)	Vinculin (19R11) Rabbit Monoclonal Antibody
描述 (Description)	Rabbit Monoclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VCL
别名 (Alternative Names)	
基因 ID (Gene ID)	7414.0
蛋白 ID (SwissProt ID)	P18206.Recombinant protein of human Vinculin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB: 1:1000
蛋白分子量 (Molecular Weight)	124kDa


研究背景 (Background)

产品名称: Vinculin (19R11) Rabbit Monoclonal Antibody
产品货号: AMRe19800



研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Vinculin in K562,C6,3T3,HeLa cell lysates using Vinculin antibody(1:1000 diluted).

注意事项 (Note)

For research use only .