

产品名称: ZNRF2 (3C15) Rabbit Monoclonal Antibody
产品货号: AMRe20293



产品概述 (Summary)

产品名称 (Production Name)	ZNRF2 (3C15) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IF-P
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNRF2
别名 (Alternative Names)	E3 ubiquitin-protein ligase ZNRF2; ZNRF2; Zinc/RING finger protein 2; RNF202; Protein Ells2;
基因 ID (Gene ID)	223082.0
蛋白 ID (SwissProt ID)	Q8NHG8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,IF-P 1:50-1:100
蛋白分子量 (Molecular Weight)	24kDa

产品名称: ZNRF2 (3C15) Rabbit Monoclonal Antibody
产品货号: AMRe20293




研究背景 (Background)

May play a role in the establishment and maintenance of neuronal transmission and plasticity via its ubiquitin ligase activity. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfer the ubiquitin to targeted substrates. May play a role in the establishment and maintenance of neuronal transmission and plasticity via its ubiquitin ligase activity. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfer the ubiquitin to targeted substrates.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ZNRF2 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .