

产品名称: SHP2 (11I13) Rabbit Monoclonal Antibody
产品货号: AMRe20370



产品概述 (Summary)

产品名称 (Production Name)	SHP2 (11I13) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% protective protein.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SHP2
别名 (Alternative Names)	Syp; Shp2; PTP1D; PTP2C; SAP-2; SHP-2; SH-PTP2; SH-PTP3; 2700084A17Rik
基因 ID (Gene ID)	19247.0
蛋白 ID (SwissProt ID)	P35235.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa; Observed MW: 68 kDa

研究背景 (Background)

Enables cell adhesion molecule binding activity; protein tyrosine phosphatase activity; and signaling receptor binding


产品名称: SHP2 (11I13) Rabbit Monoclonal Antibody
产品货号: AMRe20370



activity. Involved in negative regulation of chondrocyte differentiation; positive regulation of cytokine production; and positive regulation of ossification. Acts upstream of or within several processes, including cell surface receptor signaling pathway; myeloid cell differentiation; and regulation of hormone secretion. Predicted to be located in several cellular components, including mitochondrion; plasma membrane raft; and stress fiber. Predicted to be part of protein-containing complex. Is expressed in several structures, including alimentary system; brain; genitourinary system; hemolymphoid system gland; and liver and biliary system. Used to study several diseases, including Noonan syndrome 1; Noonan syndrome with multiple lentigines; hepatocellular adenoma; intrinsic cardiomyopathy (multiple); and juvenile myelomonocytic leukemia. Human ortholog(s) of this gene implicated in several diseases, including Noonan syndrome (multiple); Noonan syndrome with multiple lentigines 1; atrophic gastritis; juvenile myelomonocytic leukemia; and metachondromatosis. Orthologous to human PTPN11 (protein tyrosine phosphatase non-receptor type 11). [provided by Alliance of Genome Resources, Apr 2022]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa, HepG2 cells and Rat spleen, Mouse muscle tissue at 1:1000.

注意事项 (Note)

For research use only .