

**产品名称: RIP3 (Phospho S232) (16H6) Rabbit
Monoclonal Antibody**
产品货号: AMRe20391



产品概述 (Summary)

产品名称 (Production Name)	RIP3 (Phospho S232) (16H6) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RIP3 Receptor interacting protein 3; Receptor interacting serine threonine kinase 3;
别名 (Alternative Names)	Receptor interacting serine/threonine protein kinase 3; Receptor-interacting serine/threonine-protein kinase 3; RIP like protein kinase 3; RIP-3;RIP-like protein kinase 3; RIPK3;
基因 ID (Gene ID)	56532.0
蛋白 ID (SwissProt ID)	Q9QZL0.

产品应用 (Application)

**产品名称: RIP3 (Phospho S232) (16H6) Rabbit
Monoclonal Antibody**
产品货号: AMRe20391

稀释比 (Dilution Ratio)

WB 1:1000-1:2000

蛋白分子量 (Molecular Weight)


Calculated MW: 53 kDa; Observed MW: 50 kDa

研究背景 (Background)

This gene encodes a member of the POP family of proteins containing three putative transmembrane domains. This gene is expressed in cardiac and skeletal muscle and may play an important role in development of these tissues. The mouse ortholog may be involved in the regeneration of adult skeletal muscle and may act as a cell adhesion molecule in coronary vasculogenesis. Three transcript variants encoding the same protein have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Phospho-RIP3 (S232) expression in L929 cell lysate.

注意事项 (Note)

For research use only .