

产品名称: Ki-67 Rabbit Monoclonal Antibody
产品货号: AMRe20417



产品概述 (Summary)

产品名称 (Production Name)	Ki-67 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Affinity Purified

免疫原信息 (Immunogen)

基因名 (Gene Name)	MKI67
别名 (Alternative Names)	MKI67 ;Antigen KI-67
基因 ID (Gene ID)	4288.0
蛋白 ID (SwissProt ID)	P46013.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:1000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 359 kDa; Observed MW: 359 kDa

研究背景 (Background)

This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively


产品名称: Ki-67 Rabbit Monoclonal Antibody
产品货号: AMRe20417



spliced transcript variants have been described. A related pseudogene exists on chromosome X. [provided by RefSeq, Mar 2009],

研究领域 (Research Area)


图片 (Image Data)



Various whole cell lysates were separated by 4-8% SDS-PAGE, and the membrane was blotted with anti-Ki67 antibody. The HRP-


conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: HeLa Predicted band size: 359kDa

Observed band size: 359kDa




Rat colon was stained with Anti-Ki67 rabbit antibody

产品名称: Ki-67 Rabbit Monoclonal Antibody
产品货号: AMRe20417



Mouse colon was stained with Anti-Ki67 rabbit antibody



A

B

C

Immunofluorescence analysis of HEK293. Picture A: Ki67 antibody (red).Picture B: DAPI (blue). Picture C: Merge of A+B

注意事项 (Note)

For research use only .