

产品名称: Nrf2 Rabbit Monoclonal Antibody
产品货号: AMRe21024



产品概述 (Summary)

产品名称 (Production Name)	Nrf2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	NFE2L2 NRF2
别名 (Alternative Names)	Nuclear factor erythroid 2-related factor 2;NF-E2-related factor 2;NFE2-related factor 2;HEBP1;Nuclear factor, erythroid derived 2, like 2;
基因 ID (Gene ID)	4780.0
蛋白 ID (SwissProt ID)	Q16236.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000,IHC 1:1000-1:5000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:68kD;Observed MW:100kD

研究背景 (Background)


产品名称: Nrf2 Rabbit Monoclonal Antibody
产品货号: AMRe21024



Cell localization: Cytoplasmic, Nuclear. nuclear factor, erythroid 2 like 2 (NFE2L2) Homo sapiens This gene encodes a transcription factor which is a member of a small family of basic leucine zipper (bZIP) proteins. The encoded transcription factor regulates genes which contain antioxidant response elements (ARE) in their promoters; many of these genes encode proteins involved in response to injury and inflammation which includes the production of free radicals. Multiple transcript variants encoding different isoforms have been characterized for this gene. [provided by RefSeq, Sep 2015],

研究领域 (Research Area)

图片 (Image Data)



C6 whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with primary antibody(1:1000). The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .