

产品名称: Cleaved N-terminal GSDMD Rabbit Monoclonal Antibody
产品货号: AMRe21098



产品概述 (Summary)

产品名称 (Production Name)	Cleaved N-terminal GSDMD Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG, Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSDMD
别名 (Alternative Names)	
基因 ID (Gene ID)	0.0
蛋白 ID (SwissProt ID)	P57764.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000, IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:53kD; Observed MW:35kD

研究背景 (Background)


Cell localization:[Gasdermin-D]: Cytoplasm, cytosol . Inflammasome . In response to a canonical inflammasome stimulus,

产品名称: Cleaved N-terminal GSDMD Rabbit Monoclonal Antibody
产品货号: AMRe21098

such as nigericin, recruited to NLRP3 inflammasome with similar kinetics to that of uncleaved CASP1 precursor. .;
[Gasdermin-D, N-terminal]: Cell membrane ; Multi-pass membrane protein . Secreted . Released in the extracellular milieu following pyroptosis. .; [Gasdermin-D, C-terminal]: Cytoplasm, cytosol ..

研究领域 (Research Area)

图片 (Image Data)



THP-1 cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with primary antibody 1:1000. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .