

产品名称: YAP (Phospho Ser127) Rabbit Monoclonal Antibody
产品货号: AMRe21109



产品概述 (Summary)

产品名称 (Production Name)	YAP (Phospho Ser127) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Phospho
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG, Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	YAP1
别名 (Alternative Names)	YAP1; YAP65; Yorkie homolog; 65 kDa Yes-associated protein; YAP65
基因 ID (Gene ID)	10413
蛋白 ID (SwissProt ID)	P46937.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000, IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 55kD; Observed MW: 55kD

研究背景 (Background)

Cell localization: Cytoplasm, Nucleus. This gene encodes a downstream nuclear effector of the Hippo signaling pathway


产品名称: YAP (Phospho Ser127) Rabbit Monoclonal Antibody
产品货号: AMRe21109



which is involved in development, growth, repair, and homeostasis. This gene is known to play a role in the development and progression of multiple cancers as a transcriptional regulator of this signaling pathway and may function as a potential target for cancer treatment. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2013],

研究领域 (Research Area)

图片 (Image Data)



Hela cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with primary antibody 1:1000. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .