

产品名称: Paxillin (Phospho Tyr118) Rabbit Monoclonal Antibody
产品货号: AMRe21147



产品概述 (Summary)

产品名称 (Production Name)	Paxillin (Phospho Tyr118) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Phospho
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG, Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	PXN
别名 (Alternative Names)	PXN; Paxillin
基因 ID (Gene ID)	5829.0
蛋白 ID (SwissProt ID)	P49023.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000, IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 65kD; Observed MW: 65kD

研究背景 (Background)


Cell localization: Cytoplasm, Nucleus. This gene encodes a cytoskeletal protein involved in actin-membrane attachment at

产品名称: Paxillin (Phospho Tyr118) Rabbit Monoclonal Antibody
产品货号: AMRe21147

sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445). [provided by RefSeq, Aug 2011],

研究领域 (Research Area)

图片 (Image Data)



HeLa cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with primary antibody 1:1000. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .