

产品名称: HNF-3 α Rabbit Monoclonal Antibody
产品货号: AMRe21188



产品概述 (Summary)

产品名称 (Production Name)	HNF-3 α Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOXA1
别名 (Alternative Names)	FOXA1;HNF3A;TCF3A;Hepatocyte nuclear factor 3-alpha;HNF-3-alpha;HNF-3A;Forkhead box protein A1;Transcription factor 3A;TCF-3A;
基因 ID (Gene ID)	3169.0
蛋白 ID (SwissProt ID)	P55317.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000,IHC 1:200-1:1000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:49kD;Observed MW:49kD

研究背景 (Background)


产品名称: HNF-3 α Rabbit Monoclonal Antibody
产品货号: AMRe21188



Cell localization:Nucleus.This gene encodes a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver. [provided by RefSeq, Jul 2008],

研究领域 (Research Area)

图片 (Image Data)



MCF7 cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with primary antibody 1:1000. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .