

产品名称: FAM160A2 Rabbit Monoclonal Antibody
产品货号: AMRe21203



产品概述 (Summary)

产品名称 (Production Name)	FAM160A2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	FHIP1B
别名 (Alternative Names)	FHIP1B;C11orf56;FAM160A2;KIAA1759;FHF complex subunit HOOK-interacting protein 1B ;FHIP1B;FTS- and Hook-interacting protein;FHIP;
基因 ID (Gene ID)	84067.0
蛋白 ID (SwissProt ID)	Q8N612.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW;;Observed MW:107kD

研究背景 (Background)

The protein encoded by this gene is part of the FTS/Hook/FHIP (FHF) complex, which can interact with members of the


产品名称: FAM160A2 Rabbit Monoclonal Antibody
产品货号: AMRe21203



homotypic vesicular protein sorting (HOPS) complex. This interaction suggests that the encoded protein is involved in vesicle trafficking. [provided by RefSeq, Dec 2016]

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of L02 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-FAM160A2 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .