

产品名称: DPEP3 Rabbit Monoclonal Antibody
产品货号: AMRe21210



产品概述 (Summary)

产品名称 (Production Name)	DPEP3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	DPEP3;UNQ834/PRO1772
别名 (Alternative Names)	DPEP3;UNQ834/PRO1772;Dipeptidase 3;
基因 ID (Gene ID)	64180.0
蛋白 ID (SwissProt ID)	Q9H4B8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:300,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:;Observed MW:54kD

研究背景 (Background)

Cell localization:Membrane; Lipid-anchor, GPI-anchor..This gene encodes a membrane-bound glycoprotein from the family of dipeptidases involved in hydrolytic metabolism of various dipeptides, including penem and carbapenem beta-lactam


产品名称: DPEP3 Rabbit Monoclonal Antibody
产品货号: AMRe21210



antibiotics. This gene is located on chromosome 16 in a cluster with another member of this family. Alternatively spliced transcript variants that encode different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HepG2 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-DPEP3 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .