

产品名称: DHRS1 Rabbit Monoclonal Antibody
产品货号: AMRe21211



产品概述 (Summary)

产品名称 (Production Name)	DHRS1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	DHRS1;SDR19C1
别名 (Alternative Names)	DHRS1;SDR19C1;Dehydrogenase/reductase SDR family member 1 ;Short chain dehydrogenase/reductase family 19C member 1;Protein SDR19C1;
基因 ID (Gene ID)	115817.0
蛋白 ID (SwissProt ID)	Q96LJ7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000
蛋白分子量 (Molecular Weight)	Calculated MW;;Observed MW:34kD

研究背景 (Background)

Cell localization:Endoplasmic reticulum. Note=May be attached to the ER membrane by its C-terminus segment...This gene


产品名称: DHRS1 Rabbit Monoclonal Antibody
产品货号: AMRe21211



encodes a member of the short-chain dehydrogenases/reductases (SDR) family. The encoded enzyme contains a conserved catalytic domain and likely functions as an oxidoreductase. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Nov 2008],

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HepG2 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-DHRS1 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .