

产品名称: PPIL1 Rabbit Monoclonal Antibody
产品货号: AMRe21220



产品概述 (Summary)

产品名称 (Production Name)	PPIL1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPIL1;CYPL1
别名 (Alternative Names)	PPIL1;CYPL1;CGI-124;UNQ2425/PRO4984;Peptidyl-prolyl cis-trans isomerase-like 1 ;PPlase;Rotamase PPIL1;
基因 ID (Gene ID)	51645.0
蛋白 ID (SwissProt ID)	Q5TDC9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:100-1:300
蛋白分子量 (Molecular Weight)	Calculated MW;;Observed MW:19kD

研究背景 (Background)

Cell localization:Nucleus..This gene is a member of the cyclophilin family of peptidylprolyl isomerases (PPlases). The


产品名称: PPIL1 Rabbit Monoclonal Antibody
产品货号: AMRe21220



cyclophilins are a highly conserved, ubiquitous family, members of which play an important role in protein folding, immunosuppression by cyclosporin A, and infection of HIV-1 virions. Based on similarity to other PPIases, this protein could accelerate the folding of proteins and might catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HepG2 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-PPIL1 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .