

产品名称: TAF1C Rabbit Monoclonal Antibody
产品货号: AMRe21221



产品概述 (Summary)

产品名称 (Production Name)	TAF1C Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAF1C
别名 (Alternative Names)	TAF1C;TATA box-binding protein-associated factor RNA polymerase I subunit C ;RNA polymerase I-specific TBP-associated factor 110 kDa;TAFI110;TATA box-binding protein-associated factor 1C;TBP-associated factor 1C;Transcription initiation factor SL1/TIF-IB subunit C;
基因 ID (Gene ID)	9013.0
蛋白 ID (SwissProt ID)	Q15572.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW;;Observed MW:96kD

产品名称: TAF1C Rabbit Monoclonal Antibody
产品货号: AMRe21221




研究背景 (Background)

Cell localization:Nucleus, nucleolus..Initiation of transcription by RNA polymerase I requires the formation of a complex composed of the TATA-binding protein (TBP) and three TBP-associated factors (TAFs) specific for RNA polymerase I. This complex, known as SL1, binds to the core promoter of ribosomal RNA genes to position the polymerase properly and acts as a channel for regulatory signals. This gene encodes the largest SL1-specific TAF. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2011]

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HT-1080 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-TAF1C rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .