

产品名称: DNPH1 Rabbit Monoclonal Antibody
产品货号: AMRe21223



产品概述 (Summary)

产品名称 (Production Name)	DNPH1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNPH1;C6orf108;RCL
别名 (Alternative Names)	DNPH1;C6orf108;RCL;5-hydroxymethyl-dUMP N-hydrolase ;2'-deoxynucleoside 5'-phosphate N-hydrolase 1;c-Myc-responsive protein RCL;
基因 ID (Gene ID)	10591.0
蛋白 ID (SwissProt ID)	O43598.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:300
蛋白分子量 (Molecular Weight)	Calculated MW;;Observed MW:19kD

研究背景 (Background)

Cell localization:Cytoplasm. Nucleus..This gene was identified on the basis of its stimulation by c-Myc protein. The latter is a


产品名称: DNPH1 Rabbit Monoclonal Antibody
产品货号: AMRe21223



transcription factor that participates in the regulation of cell proliferation, differentiation, and apoptosis. The exact function of this gene is not known but studies in rat suggest a role in cellular proliferation and c-Myc-mediated transformation. Two alternative transcripts encoding different proteins have been described. [provided by RefSeq, Jul 2008],

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of K-562 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-DNPH1 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .