

产品名称: POC5 Rabbit Monoclonal Antibody
产品货号: AMRe21224



产品概述 (Summary)

产品名称 (Production Name)	POC5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	POC5;C5orf37
别名 (Alternative Names)	POC5;C5orf37;Centrosomal protein POC5 ;Protein of centriole 5;hPOC5;
基因 ID (Gene ID)	134359.0
蛋白 ID (SwissProt ID)	Q8NA72.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:300
蛋白分子量 (Molecular Weight)	Calculated MW;;Observed MW:63kD

研究背景 (Background)

Cell localization: Cytoplasm, cytoskeleton, microtubule organizing center, centrosome. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole. Note=Localized to the distal portion of centrioles..Predicted to be involved in eye


产品名称: POC5 Rabbit Monoclonal Antibody
产品货号: AMRe21224



photoreceptor cell development. Located in centrosome; cytosol; and nucleoplasm. [provided by Alliance of Genome Resources, Feb 2025]

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HEK-293 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-POC5 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .