

产品名称: LYPD3 Rabbit Monoclonal Antibody
产品货号: AMRe21232



产品概述 (Summary)

产品名称 (Production Name)	LYPD3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	LYPD3;C4.4A
别名 (Alternative Names)	LYPD3;C4.4A;UNQ491/PRO1007;Ly6/PLAUR domain-containing protein 3 ;GPI-anchored metastasis-associated protein C4.4A homolog;Matrigel-induced gene C4 protein;MIG-C4;
基因 ID (Gene ID)	27076.0
蛋白 ID (SwissProt ID)	O95274.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:300,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW::Observed MW:38kD

研究背景 (Background)


产品名称: LYPD3 Rabbit Monoclonal Antibody
产品货号: AMRe21232



Cell localization: Cell membrane; Lipid-anchor, GPI-anchor..Predicted to enable laminin binding activity. Involved in negative regulation of smooth muscle cell apoptotic process. Located in extracellular space. [provided by Alliance of Genome Resources, Feb 2025]

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HepG2 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-LYPD3 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .