

产品名称: Prostein Rabbit Monoclonal antibody
产品货号: AMRe21233



产品概述 (Summary)

产品名称 (Production Name)	Prostein Rabbit Monoclonal antibody
描述 (Description)	Rabbit Monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC45A3
别名 (Alternative Names)	SLC45A3;PCANAP6;PRST;Solute carrier family 45 member 3 ;Prostate cancer-associated protein 6;Prostein;
基因 ID (Gene ID)	85414.0
蛋白 ID (SwissProt ID)	Q96JT2.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-300
蛋白分子量 (Molecular Weight)	Calculated MW;;Observed MW:61kD

研究背景 (Background)

Cell localization:Membrane; Multi-pass membrane protein..Predicted to enable sucrose:proton symporter activity. Predicted


产品名称: Prostein Rabbit Monoclonal antibody
产品货号: AMRe21233



to be involved in positive regulation of small molecule metabolic process; regulation of oligodendrocyte differentiation; and sucrose transport. Predicted to be located in plasma membrane. Predicted to be active in membrane. [provided by Alliance of Genome Resources, Apr 2022]

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human Prostate tissue. 1,primary Antibody was diluted at 1:200(4°C,overnight). 2, EDTA pH 9.0 was used for antibody retrieval(>98°C,20min). 3,Secondary antibody was diluted at 1:200(room temperature, 30min).

注意事项 (Note)

For research use only .