

产品名称: RBMY1A1 Rabbit Monoclonal Antibody
产品货号: AMRe21239



产品概述 (Summary)

产品名称 (Production Name)	RBMY1A1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBMY1A1;RBM1
别名 (Alternative Names)	RBMY1A1;RBM1;RBM2;YRRM1;YRRM2;RNA-binding motif protein, Y chromosome, family 1 member A1 ;RNA-binding motif protein 1;RNA-binding motif protein 2;Y chromosome RNA recognition motif 1;hRBMY;
基因 ID (Gene ID)	5940.0
蛋白 ID (SwissProt ID)	P0DJD3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:300
蛋白分子量 (Molecular Weight)	Calculated MW:;Observed MW:55kD

研究背景 (Background)


产品名称: RBMY1A1 Rabbit Monoclonal Antibody
产品货号: AMRe21239



Cell localization:Nucleus..This gene encodes a protein containing an RNA-binding motif in the N-terminus and four SRGY (serine, arginine, glycine, tyrosine) boxes in the C-terminus. This protein is thought to function as a splicing regulator during spermatogenesis. Multiple closely related paralogs of this gene are found in a gene cluster in the AZFb azoospermia factor region of chromosome Y. Most of these related copies are thought to be pseudogenes, though several likely encode functional proteins. [provided by RefSeq, Mar 2016]

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of COLO 320 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-RBMY1A1 rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .