

产品名称: TAGAP Rabbit Monoclonal Antibody
产品货号: AMRe21240



产品概述 (Summary)

产品名称 (Production Name)	TAGAP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAGAP;TAGAP1;FKSG15
别名 (Alternative Names)	TAGAP;TAGAP1;FKSG15;T-cell activation Rho GTPase-activating protein ;T-cell activation GTPase-activating protein;
基因 ID (Gene ID)	117289.0
蛋白 ID (SwissProt ID)	Q8N103.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,FC 1:100-1:300,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW;;Observed MW:80kD

研究背景 (Background)

This gene encodes a member of the Rho GTPase-activator protein superfamily. The encoded protein may function as a Rho


产品名称: TAGAP Rabbit Monoclonal Antibody
产品货号: AMRe21240



GTPase-activating protein. Alterations in this gene may be associated with several diseases, including rheumatoid arthritis, celiac disease, and multiple sclerosis. Alternate splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2013],

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HEK-293 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-TAGAP rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .