

产品名称: CPVL Rabbit Monoclonal Antibody
产品货号: AMRe21244



产品概述 (Summary)

产品名称 (Production Name)	CPVL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	CPVL;VLP;PSEC0124;UNQ197/PRO223 CPVL;VLP;PSEC0124;UNQ197/PRO223;Probable serine carboxypeptidase
别名 (Alternative Names)	CPVL ;Carboxypeptidase, vitellogenic-like;Vitellogenic carboxypeptidase-like protein;VCP-like protein;hVLP;
基因 ID (Gene ID)	54504.0
蛋白 ID (SwissProt ID)	Q9H3G5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:300
蛋白分子量 (Molecular Weight)	Calculated MW::Observed MW:52kD

研究背景 (Background)


产品名称: CPVL Rabbit Monoclonal Antibody
产品货号: AMRe21244



The protein encoded by this gene is a carboxypeptidase and bears strong sequence similarity to serine carboxypeptidases. Carboxypeptidases are a large class of proteases that act to cleave a single amino acid from the carboxy termini of proteins or peptides. The exact function of this protein, however, has not been determined. [provided by RefSeq, Jan 2017]

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HEK-293 whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-CPVL rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody.

注意事项 (Note)

For research use only .