

产品名称: β Tubulin Rabbit Monoclonal Antibody
产品货号: AMRe21344



产品概述 (Summary)

产品名称 (Production Name)	β Tubulin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	TUBB
别名 (Alternative Names)	Tubulin beta chain (Tubulin beta-5 chain)
基因 ID (Gene ID)	203068
蛋白 ID (SwissProt ID)	P07437.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000,IHC 1:2000-1:10000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:50kD;Observed MW:50kD

研究背景 (Background)

Cell localization:Cytoplasm.This gene encodes a beta tubulin protein. This protein forms a dimer with alpha tubulin and acts


产品名称: β Tubulin Rabbit Monoclonal Antibody
产品货号: AMRe21344



as a structural component of microtubules. Mutations in this gene cause cortical dysplasia, complex, with other brain malformations 6. Alternative splicing results in multiple splice variants. There are multiple pseudogenes for this gene on chromosomes 1, 6, 7, 8, 9, and 13. [provided by RefSeq, Jun 2014],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Hela cells, using β Tubulin Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .