

产品名称: ATP5A Rabbit Monoclonal Antibody
产品货号: AMRe21352



产品概述 (Summary)

产品名称 (Production Name)	ATP5A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATP5A1
别名 (Alternative Names)	ATP5A1;ATP5A;ATP5AL2;ATPM;ATP synthase subunit alpha; mitochondrial
基因 ID (Gene ID)	498
蛋白 ID (SwissProt ID)	P25705.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000,IHC 1:200-1:1000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:60kD;Observed MW:55kD

研究背景 (Background)

Cell localization:Mitochondrion.This gene encodes a subunit of mitochondrial ATP synthase. Mitochondrial ATP synthase


产品名称: ATP5A Rabbit Monoclonal Antibody
产品货号: AMRe21352



catalyzes ATP synthesis, using an electrochemical gradient of protons across the inner membrane during oxidative phosphorylation. ATP synthase is composed of two linked multi-subunit complexes: the soluble catalytic core, F1, and the membrane-spanning component, Fo, comprising the proton channel. The catalytic portion of mitochondrial ATP synthase consists of 5 different subunits (alpha, beta, gamma, delta, and epsilon) assembled with a stoichiometry of 3 alpha, 3 beta, and a single representative of the other 3. The proton channel consists of three main subunits (a, b, c). This gene encodes the alpha subunit of the catalytic core. Alternatively spliced transcript variants encoding the different isoforms have been identified. Pseudogenes of thi

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from A549

cells, using ATP5A Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .