

产品名称: VWF Rabbit Monoclonal Antibody
产品货号: AMRe21354



产品概述 (Summary)

产品名称 (Production Name)	VWF Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	VWF
别名 (Alternative Names)	F8VWF
基因 ID (Gene ID)	7450.0
蛋白 ID (SwissProt ID)	P04275.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000,IHC 1:1000-1:4000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:309kD;Observed MW:280kD

研究背景 (Background)

Cell localization:Secreted.This gene encodes a glycoprotein involved in hemostasis. The encoded preproprotein is


产品名称: VWF Rabbit Monoclonal Antibody
产品货号: AMRe21354



proteolytically processed following assembly into large multimeric complexes. These complexes function in the adhesion of platelets to sites of vascular injury and the transport of various proteins in the blood. Mutations in this gene result in von Willebrand disease, an inherited bleeding disorder. An unprocessed pseudogene has been found on chromosome 22.
[provided by RefSeq, Oct 2015],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HepG2

cells, using VWF Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .