

产品名称: Ly6g Rabbit Monoclonal Antibody
产品货号: AMRe21366



产品概述 (Summary)

产品名称 (Production Name)	Ly6g Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	Ly6g
别名 (Alternative Names)	Gr-1;Gr1;Ly-6G.1;Ly6G;Lymphocyte antigen 6 complex locus G;Lymphocyte antigen 6G
基因 ID (Gene ID)	546644.0
蛋白 ID (SwissProt ID)	P35461.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:1000-1:4000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:14kD;Observed MW:14kD

研究背景 (Background)


产品名称: Ly6g Rabbit Monoclonal Antibody
产品货号: AMRe21366



Cell localization:Membrane.Predicted to enable acetylcholine receptor binding activity and acetylcholine receptor inhibitor activity. Predicted to be involved in acetylcholine receptor signaling pathway. Located in external side of plasma membrane. Is expressed in large intestine; small intestine; and spleen. [provided by Alliance of Genome Resources, Nov 2024]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Mouse kidney

cells, using Ly6g Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .