

产品名称: Bcl-2 Rabbit Monoclonal Antibody  
产品货号: AMRe21374



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | Bcl-2 Rabbit Monoclonal Antibody       |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA,IP                 |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                        |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | IgG,Kappa                                                                                |
| 克隆 (Clonality)      | Monoclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein                           |
| 纯化方式 (Purification) | Protein A                                                                                |

## 免疫原信息 (Immunogen)

|                        |                                |
|------------------------|--------------------------------|
| 基因名 (Gene Name)        | BCL2                           |
| 别名 (Alternative Names) | BCL2;Apoptosis regulator Bcl-2 |
| 基因 ID (Gene ID)        | 596                            |
| 蛋白 ID (SwissProt ID)   | P10415.                        |

## 产品应用 (Application)

|                          |                                                                                           |
|--------------------------|-------------------------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:2000-1:10000,IHC 1:500-1:2000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200 |
| 蛋白分子量 (Molecular Weight) | Calculated MW:26kD;Observed MW:26kD                                                       |

## 研究背景 (Background)

Cell localization:Cytoplasm.BCL2, apoptosis regulator(BCL2) Homo sapiens This gene encodes an integral outer


**产品名称: Bcl-2 Rabbit Monoclonal Antibody**  
**产品货号: AMRe21374**



mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016],

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of lysates from MCF7 cells, using Bcl-2 Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

## 注意事项 (Note)

For research use only .