

产品名称: c-Jun Rabbit Monoclonal Antibody
产品货号: AMRe21376



产品概述 (Summary)

产品名称 (Production Name)	c-Jun Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	JUN
别名 (Alternative Names)	Transcription factor AP-1 (Activator protein 1;AP1;Proto-oncogene c-Jun;V-jun avian sarcoma virus 17 oncogene homolog;p39)
基因 ID (Gene ID)	3725.0
蛋白 ID (SwissProt ID)	P05412.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000,IHC 1:200-1:1000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:36kD;Observed MW:43kD

研究背景 (Background)


产品名称: c-Jun Rabbit Monoclonal Antibody
产品货号: AMRe21376



Cell localization:Nucleus.function:Transcription factor that recognizes and binds toThe enhancer heptamer motif 5'-TGA[CG]TCA-3',PTM:Phosphorylation enhancesThe transcriptional activity. Phosphorylated by PRKDC.,similarity:Belongs toThe bZIP family.,similarity:Belongs toThe bZIP family. Jun subfamily.,similarity:Contains 1 bZIP domain.,subunit:Heterodimer with either FOS or BATF3. Interacts with HIVP3 (By similarity). Interacts with SMAD3/SMAD4 heterodimers. Interacts with MYBBP1A, SPIB and TCF20. Interacts with COPS5; indirectly leading to its phosphorylation. Interacts with DSIPI; this interaction inhibitsThe binding of active AP1 to its target DNA.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from PC-12

cells, using c-Jun Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .