

产品名称: Cardiac Troponin I Rabbit Monoclonal Antibody
产品货号: AMRe21381



产品概述 (Summary)

产品名称 (Production Name)	Cardiac Troponin I Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNNI3
别名 (Alternative Names)	TNNI3;TNNC1;Troponin I; cardiac muscle;Cardiac troponin I
基因 ID (Gene ID)	7137
蛋白 ID (SwissProt ID)	P19429.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000,IHC 1:200-1:1000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:24kD;Observed MW:24kD

研究背景 (Background)


Cell localization:Cytoplasm.Troponin I (TnI), along with troponin T (TnT) and troponin C (TnC), is one of 3 subunits that

产品名称: Cardiac Troponin I Rabbit Monoclonal Antibody
产品货号: AMRe21381

formThe troponin complex ofThe thin filaments of striated muscle. TnI isThe inhibitory subunit; blocking actin-myosin interactions and thereby mediating striated muscle relaxation.The TnI subfamily contains three genes: TnI-skeletal-fast-twitch, TnI-skeletal-slow-twitch, and TnI-cardiac. This gene encodesThe TnI-cardiac protein and is exclusively expressed in cardiac muscle tissues. Mutations in this gene cause familial hypertrophic cardiomyopathy type 7 (CMH7) and familial restrictive cardiomyopathy (RCM). [provided by RefSeq, Jul 2008],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Rat heart

cells, using Cardiac Troponin I Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .